

COMMODITY SPECIFICATION

Product	Airway Management Simulator	Model	BT-CSIE
Manufacture	BT INC	Q'ty	1 SET



A. FEATURES

- Objective intubation practice with an adult torso model
- Ability to do realistic practice by using skin-like feeling Human torso model
- Head tilt / Chin lift, Sniffing position
- Jaw Thrust maneuver (Displayed through the monitor)
- Ventilation by using Bag-Valve-Mask
- Insert the oral and nasal airway can be applied
- Intubation tube (ETT), laryngeal mask (LMA/l-gel), combi tube (Combi-Tube) can be applied
- Intubation by using laryngoscope
- Observable chest rising during the ventilation and check with stethoscope
- Alarm sounds due to the excessive use of laryngoscope
- Detects Head tilt-chin lift method, Jaw Thrust, sniffing position and display these through a tablet
- Observable intubation depth (deep/moderate/shallow) through tablet
- Observable esophageal intubation on the tablet PC
- Detects tidal volume and stomach expansion and display the state on the Tablet PC
- User can check up the training record through the auto-ID items and check-list
- Check-list & result-point is all configurable.
- All the results can be saved.
- Wireless Blue-tooth connection between display device-tablet PC

COMMODITY SPECIFICATION

- Power : AC power or battery

B. SPECIFICATION & CONSIST OF

- * Adult Torso Model – 1EA (Controller(installed) / Blue-tooth device(installed))
- * Tablet PC – 1EA (Software include)
- * Lubricant 1EA
- * AC Power Adaptor..... 1EA
- * Battery Holder..... 1EA
- * Storage Case..... 1EA
- * User Guide 1EA

C. WEIGHT & SIZE

- Adult Torso Manikin : 8.5kg (700×300×220mm)
- Display device : 0.6kg (263×181×10mm)

COMMODITY SPECIFICATION

--