

# COMMODITY SPECIFICATION

|                    |                                       |              |                |
|--------------------|---------------------------------------|--------------|----------------|
| <b>Product</b>     | <b>Enema Administration Simulator</b> | <b>Model</b> | <b>BT-CSEN</b> |
| <b>Manufacture</b> | <b>BT INC</b>                         | <b>Q'ty</b>  | <b>1 SET</b>   |



## A. FEATURES

- Life-like practical training with BT Product
- Make the Sim's Position, Knee-chest Position
- Any water-soluble Jelly being used at clinical situation can be applied
- Any Glycerin being used at clinical situation can be applied
- Training with various size of catheter is available
- Non-return valve : Prevent leaking after taking catheter out
- Artificial feces can be inserted / removed

## B. CONSIST OF

- Adult Buttocks model : 1set
- Drain tube : 1EA
- Storage bag : 1EA

## C. WEIGHT & SIZE

- Adult Buttocks model : 3.5kg (330×430×340mm)
- Storage bag : 2.4kg (440×330×390mm)
- Total : Approximately 6kg