


COMMODITY SPECIFICATION

| | | | |
|--------------------|--|--------------|----------------|
| Product | Female Catheterization Training Model | Model | BT-CSCF |
| Manufacture | BT INC | Q'ty | 1 EA |



A. FEATURES

- Life-like practical training
- Make position similar to the real patient
- Any water-soluble Jelly being used at clinical situation can be applied
- Training with various size of catheter is available
- Realistic practice with attachable extended legs
- Non-return valve : Prevent leaking after taking catheter out.
- Checking fluid after inserting a catheter into bladder correctly
- Alarm when a catheter inserted into a vagina ; to operate, a battery is needed
- Genital organs and buttocks are attachable : easy to maintain

B. CONSIST OF

- Adult Buttocks model : 1sets
- Extended legs : 1sets
- Simulated urine bag : 1EA
- Storage case : 1EA

COMMODITY SPECIFICATION

C. WEIGHT & SIZE

- Adult buttocks model : 3.5kg (430×280×200mm)
- Extended legs (Left/Right) : Each 1kg (320×150×150mm)
- Storage case : 2.4kg (440×330×390mm)
- Total : Approximately 8kg