

COMMODITY SPECIFICATION

Product	Airway Management Training Model	Model	BT-CSIE-T
Manufacture	BT INC	Q'ty	1 SET



A. FEATURES

- Objective intubation practice with an adult torso model
- Ability to do realistic practice by using skin-like feeling Human torso model
- Head tilt / Chin lift, Sniffing position
- Jaw Thrust maneuver
- Ventilation by using Bag-Valve-Mask
- Insert the oral and nasal airway can be applied
- Intubation tube (ETT), laryngeal mask (LMA/l-gel), combi tube (Combi-Tube) can be applied
- Intubation by using laryngoscope
- Observable chest rising during the ventilation and check with stethoscope
- Observable gastric inflation during ventilation while the ETT is inserted into a stomach

COMMODITY SPECIFICATION

B. SPECIFICATION & CONSIST OF

- * Adult Torso Model – 1EA
- * Storage Case – 1EA
- * User Guide – 1EA

C. WEIGHT & SIZE

- Adult Torso Manikin : 8.5kg (700×300×240mm)