

COMMODITY SPECIFICATION

Product	Airway Management Simulator	Model	BT-CSIE
Manufacture	BT INC	Q'ty	1 SET



A. FEATURES

- Objective intubation practice with an adult torso model
- Realistic training with life-like skin
- Head tilt / Chin lift, Sniffing position (Displayed through the monitor)
- Jaw Thrust maneuver (Displayed through the monitor)
- Ventilation with a Bag-Valve-Mask
- Insert the oral airway (OPA) and nasal airway(NPA) to the manikin
- Intubation tube (ETT), laryngeal mask (LMA/I-gel), combi tube (Combi-Tube) can be inserted
- Intubation by using laryngoscope
- Observable chest rising and check the sound with stethoscope during the ventilation
- Alarm sounds for the front teeth due to the excessive use of laryngoscope
- Detects Head tilt-chin lift method, Jaw Thrust, sniffing position and display these through a tablet
- Observable intubation depth (deep/moderate/shallow) through tablet
- Observable esophageal intubation on the tablet PC
- Detects tidal volume, gastric inflation and display the state on the Tablet PC
- User can check up the training record through the auto-check list and manual check list
- Adjust each score for all the check-list items
- All the results can be saved.
- Wireless Blue-tooth connection between display device-tablet PC
- Power : AC power or battery; can be used wirelessly

COMMODITY SPECIFICATION

B. SPECIFICATION & CONSIST OF

- * Adult Torso Model – 1EA (Controller(installed) / Blue-tooth device(installed))
- * Tablet PC – 1EA (Software include)
- * Lubricant 1EA
- * AC Power Adaptor..... 1EA
- * Battery Holder..... 1EA
- * Storage Case..... 1EA
- * User Guide 1EA

C. WEIGHT & SIZE

- Adult Torso Manikin : 8.5kg (700×300×240mm)