

COMMODITY SPECIFICATION

| | | | |
|--------------------|--|--------------|----------------|
| Product | Tube Feeding & Tracheostomy care & Suction Training Model | Model | BT-CSIS |
| Manufacture | BT INC | Q'ty | 1 SET |



A. FEATURES

- Realistic training with anatomically correct structures.
- Oral/Nasal feeding tube injection.
- Check tube placement by auscultation with stethoscope.
- Practicing with a real liquid food
- Easy to maintain ; Detachable stomach and trachea.
- Positioning with a stand by adjusting the angle. (0°, 21~45°, 90°)
- Checking tube placement through transparent window
- Training tube feeding, PEG care, Tracheostomy care, suction, gastrostomy tube care

COMMODITY SPECIFICATION

B. SPECIFICATION & CONSIST OF

- Adult upper torso : 1set
- Drain Tube : 1set
- Storage case : 1EA

C. WEIGHT & SIZE

- Adult upper torso : 7.8kg (310×700×238mm)
- Storage case : 7kg (380×820×270mm)
- Total : Approximately 15kg