

COMMODITY SPECIFICATION

| | | | |
|--------------------|---|--------------|------------------|
| Product | Self-training CPR Model(Sherpa Plus) | Model | BT-SEEM-P |
| Manufacture | BT INC. | Q'ty | 1 SET |



A. FEATURES

- Chest Compression, Airway, Ventilation skill practice needed for basic life support available.
- Total number of chest compression and number of good chest compression are displayed.
- Instant compression depth displayed.
- Chest compression depth and rate can be checked with LED lamp.
- Detect Incomplete recoil and over-compression.
- Sound beeps once user make the proper compression depth and speed.
- Sound and LED Lamp can be turned on/off for an advanced training.
- Checking chest rising visually once user gives the proper artificial respiration.
- By attaching offering magnetic stickers to AED trainer pads, AED trainer pads can be used semi-permanently.
- Hygienic training condition with self-repairable face-skin and lung.
- Storage case can be used as a mattress pad.
- Operated with Batteries and AC adapter.
- Can be upgraded to Sherpa X with a purchase of the upgrade kit when customers want.

COMMODITY SPECIFICATION

B. COMPONENTS

- Adult Torso : 1EA
- Storage case : 1EA
- AC Power Adaptor : 1EA
- User Manual : 1EA

C. WEIGHT & SIZE

- Adult Torso : 4.4kg (650×335×215mm)
- Storage case : 1.5kg (670×350×230mm)
- Total : Approximately 6.2kg