



COMMODITY SPECIFICATION

Product	SMART CPR Training Model (Sherpa X)	Model	BT-SHERPA-XM
Manufacture	BT INC / South Korea	Q'ty	1 SET



A. FEATURES

1. SMART CPR training model(Sherpa X), Training & Evaluation up to 1:30

- Connect up to 6 manikins to 1 tablet through Bluetooth
- Tablet and laptop are connected with WiFi, Possible to train & evaluate up to 1:30
- Easy to use and user friendly GUI software
- CPR training and evaluation through its own S/W
- Saving and printing the training result through its own S/W
- Evaluation result can be saved as JPEG or Excel format
- Change layout to 1:12, 1:18, 1:24, 1:30 according to the system's setting
- Set the guideline
- Group setting and typing the user name
- On/Off the LED display panel of the manikin
- Make a statistics and analyze the results

2. A case set for Multi-system

- Charge & Store tablets
- Wifi router is included
- Make sure the safety with a lock

COMMODITY SPECIFICATION

3. CPR Training Model ; Sherpa X

- **Wireless connection with a tablet or a smartphone through Bluetooth**
 - Easy to use and user friendly GUI software
 - Connect up to 6 manikins to 1 tablet (tablet is optional)
 - Compression depth & rate and hand position, ventilation volume can be checked through its own S/W
 - Saving and printing the training result through its own S/W
 - Scoring the training result and rank the result in order of high scores
 - Select the mode between Training mode & Evaluation mode

- **LED Display embedded in the manikin**
 - Real-time feedback of chest compression depth and rate with LED lamp
 - Display total number of compression and the number of good compression
 - Display the instant depth of compression
 - Detect incomplete recoil
 - Detect the good chest compression when depth and rate are good and give sound feedback (Sound On/Off)
 - Display the feedback about ventilation such as volume and expiratory time
 - Display the total number of ventilation and the number of good ventilation
 - Display the ventilation volume

- **Manikin Body**
 - Observable chest rising during ventilation
 - By attaching offering magnetic stickers to AED trainer pads, AED trainer pads can be used semi-permanently
 - Easy to replace the face skin and lung ; sanitary training environment
 - The storage bag can be used as a training mattress

COMMODITY SPECIFICATION

B. COMPONENTS

- Laptop with software : 1set
- Case set for Multi-system : 1set (Wifi router, Multi-charger, Extension cord, Cabinet)
- SMART CPR training model : Optional (Depends on the Multi-system's composition)
- 10.1" tablet : Optional (Depends on the Multi-system's composition)

C. WEIGHT & SIZE

- Laptop : Depends on each laptop's specification
- Case set for Multi-system : 20kg (500×600×500mm)