

COMMODITY SPECIFICATION

Product	SMART CPR Training Model (Sherpa X)	Model	BT-SEEM2
Manufacture	BT INC.	Q'ty	1 SET



A. FEATURES

- Chest Compression, Airway, Ventilation skill practice needed for basic life support available.
- Total number of chest compression and number of good chest compression are displayed.
- Instant compression depth displayed.
- Chest compression depth and rate can be checked with LED lamp.
- Detect Incomplete recoil and over-compression.
- Sound beeps once user make the proper compression depth and speed.
- Sound and LED Lamp can be turned on/off for an advanced training.
- Checking chest rising visually once user gives the proper artificial respiration.
- By attaching offering magnetic stickers to AED trainer pads, AED trainer pads can be used semi-permanently.
- Hygienic training condition with self-repairable face-skin and lung.
- Storage case can be used as a mattress pad.
- Operated with Batteries and AC adapter.
- Detailed information about chest compression(depth, location, rate, number, recoil, hands-off time and etc.) and ventilation(volume, rate, number and average etc.) can be displayed through BT's S/W.
- Compression location, ventilation can be checked through BT's S/W.
- Training and evaluation with BT's S/W.
- Available to check evaluation result and print it out through BT's S/W.

COMMODITY SPECIFICATION

- Evaluation up to 6 manikins through Bluetooth connection.
- Able to evaluation a group and make a rank by using a training result.
- Training result can be saved in PNG format and exported to Excel format.
- Guideline can be adjusted through S/W.

B. COMPONENTS

- Adult Torso : 1EA
- Storage case : 1EA
- AC Power Adaptor
- User Manual : 1EA

C. WEIGHT & SIZE

- Adult Torso : 4.4kg (650×335×215mm)
- Storage case : 1.5kg (670×350×230mm)
- Total : Approximately 6.2kg